

Tissue Cassette

Labsim tissue cassette: Printing information is high temperature resistant, corrosion resistant and compatible with the mainstream tissue cassette printers in the market. It has the characteristics of clear and complete writing, long-lasting information storage and non-fading, especially suitable for large-scale pathological diagnosis. Support information scanning anytime and anywhere, which facilitates quick identification of paraffin block and improves work efficiency.

No.	Item No.	Product Name	Description	Picture
1	Q.Q2.000144	Sponge for tissue cassette (white)	box (10 bags × 1000 pcs)	
	Q.Q2.000145	Sponge for tissue cassette (blue)		
2	W.W5.000471	Thermal tissue cassette (white)	Hole size: 4.7 × 0.75 mm Outer dimension with lid: 41.3(L) × 28.2(W) × 6(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 28.7(L) × 25.7(W) × 3.7(H) mm <hr/> 75 pcs / tray, 20 trays / box 80 trays / carton, (w / o lid, Cover Lid will be packed separately)	
	W.W5.000472	Thermal tissue cassette (pink)		
	W.W5.000474	Thermal tissue cassette (yellow)		
	W.W5.000475	Thermal tissue cassette (blue)		
	W.W5.000476	Thermal tissue cassette (green)		
	W.W5.000477	Thermal tissue cassette (purple)		
3	Q.Q2.000119	Laser tissue cassette (white)	Hole size: 5 × 1 mm Outer dimension with lid: 43(L) × 28.2(W) × 7(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 39.9(L) × 26.2(W) × 5.1(H) mm <hr/> 75 pcs / tray, 30 trays / carton, (with lid)	
	Q.Q2.000120	Laser tissue cassette (pink)		
	Q.Q2.000121	Laser tissue cassette (green)		
	Q.Q2.000122	Laser tissue cassette (yellow)		
	Q.Q2.000123	Laser tissue cassette (purple)		
	Q.Q2.000124	Laser tissue cassette (blue)		
4	Q.Q2.000126	Laser tissue cassette (white)	Hole size: 5 × 1 mm Outer dimension with lid: 43(L) × 28.2(W) × 7(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 39.9(L) × 26.2(W) × 5.1(H) mm <hr/> 90 pcs / tray 20 trays / carton (w / o lid, cover lid will be pack separately)	
	Q.Q2.000127	Laser tissue cassette (pink)		
	Q.Q2.000128	Laser tissue cassette (yellow)		
	Q.Q2.000129	Laser tissue cassette (green)		
	Q.Q2.000130	Laser tissue cassette (purple)		
	Q.Q2.000131	Laser tissue cassette (blue)		

No.	Item No.	Product Name	Description	Picture
5	Q.Q2.000153	Laser tissue cassette with cover (white)	Hole size: 0.5 × 0.5 mm Outer dimension with lid: 41.3(L) × 28.2(W) × 6(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 28.7(L) × 25.7(W) × 3.7(H) mm <hr/> 85 pcs / tray 4 trays / boxes 136trays / carton (with lid)	
	Q.Q2.000154	Laser tissue cassette with cover (pink)		
	Q.Q2.000155	Laser tissue cassette with cover (yellow)		
	Q.Q2.000156	Laser tissue cassette with cover (purple)		
	Q.Q2.000157	Laser tissue cassette with cover (green)		
	Q.Q2.000158	Laser tissue cassette with cover (blue)		
6	Q.Q2.000167	Laser tissue cassette separate with cover (white)	Hole size: 0.5 × 0.5 mm Outer dimension with lid: 41.3(L) × 28.2(W) × 6(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 28.7(L) × 25.7(W) × 3.7(H) mm <hr/> 85 pcs / tray 4 trays / boxes 136trays / carton (w / o lid, cover lid will be pack separately)	
	Q.Q2.000168	Laser tissue cassette separate with cover (pink)		
	Q.Q2.000169	Laser tissue cassette separate with cover (yellow)		
	Q.Q2.000170	Laser tissue cassette separate with cover (purple)		
	Q.Q2.000171	Laser tissue cassette separate with cover (green)		
	Q.Q2.000172	Laser tissue cassette separate with cover (blue)		
7	Q.Q2.000181	Laser tissue cassette with cover (white)	Hole size: 0.26 × 0.26 mm Outer dimension with lid: 41.3(L) × 28.2(W) × 6(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 28.7(L) × 25.7(W) × 3.7(H) mm <hr/> 75 pcs / tray 72 trays / carton (with lid)	
	Q.Q2.000182	Laser tissue cassette with cover (pink)		
	Q.Q2.000183	Laser tissue cassette with cover (yellow)		
	Q.Q2.000184	Laser tissue cassette with cover (purple)		
	Q.Q2.000185	Laser tissue cassette with cover (green)		
	Q.Q2.000186	Laser tissue cassette with cover (blue)		
8	Q.Q2.000195	Laser tissue cassette separate with cover (white)	Hole size: 0.26 × 0.26 mm Outer dimension with lid: 41.3(L) × 28.2(W) × 6(H) mm Outer dimension without lid: 41.3(L) × 28.2(W) × 6(H) mm Inner dimension: 28.7(L) × 25.7(W) × 3.7(H) mm <hr/> 75 pcs / tray 72 trays / carton (w / o lid, cover lid will be pack separately)	 
	Q.Q2.000196	Laser tissue cassette separate with cover (pink)		
	Q.Q2.000197	Laser tissue cassette separate with cover (yellow)		
	Q.Q2.000198	Laser tissue cassette separate with cover (purple)		
	Q.Q2.000199	Laser tissue cassette separate with cover (green)		
	Q.Q2.000200	Laser tissue cassette separate with cover (blue)		

Note: MOQ is 5 cartons