



# G364

## Slide Laser Printer

Slide laser printer is specially designed for standard slide printing through laser technology.

Easy to operate, printed content is in high quality, corrosion resistant and suitable for long-term preservation.



### Clear Anti-corrosion

The printed text is clear, resistant to xylene, alcohol and other solvents corrosion



### Highly com- patible slides

Can print a variety of paint size slides, perfect for different color slides



### New UV laser

The laser spot is extremely small and precision is higher, and can print extremely high-resolution QR codes and related information clearly



### Integrated design

Independent integrated design, stand-alone printing without external computer

## Technical parameters

- Technology: UV laser printing
  - Data transmission : USB, Ethernet, serial ports
  - Loading capacity: 100 pcs
  - Output collection capacity:
    - Slide storage box: 10 pcs (Print and Use)
    - Big storage box: 35 / 50 / 80 pcs (Batch Printing)
  - Printing accuracy: 2500 dpi
  - Print speed: ≤ 3 s / piece
  - Size: 200 (W) \* 437 (D) \* 445 (H) mm
- Built-in barcode scanner and external barcode scanner
  - Non-contact laser technology without consumables
  - 10.1-inch full-color LCD touch screen
  - Clearly print characters, letters, symbols, barcodes and other information
  - The imprints are corrosion resistant exposed to xylene, alcohol and decalcifying fluid, enabling long term preservation
  - Easy to operate with professional slide loading / output device
  - Built-in operating system without external connection
  - System: LIS/HIS/PIS/PACS and other hospital commonly used systems are compatible
  - Composite air treatment function



---

**Contact Us**  
**Zenith Biomedical**

📍 1265 Barry Ave, PH2, Los Angeles, CA 90025

☎ 612-433-7719

✉ [info@zenithbiomedical.com](mailto:info@zenithbiomedical.com)