

— Make Pathology Easier

TB406I Tissue Cassette Laser Printer

Designed for standard tissue cassette printing, it uses laser technology to print tissue cassettes. The instrument is easy to operate, and the printed content is high quality and corrosion-resistant, suitable for long-term storage.



Small • Intelligent • Convenient





High-quality printing

The printed content is clear and can withstand water and xylene. It can be soaked in alcohol and other solutions and is suitable for long-term storage

The template can be edited freely

Free drag and drop editing, Arbitrary printing in the tissue cassette printing area, support for multiple templates



Barcode scanning module

Equipped with a barcode scanning module, the barcode reader can be rotated to a comfortable angle, and the software controls the switch of the barcode scanning module



Exquisite and compact

The instrument is small in size, with a width of $\leq 11\text{cm}$, which greatly saves space



Data transmission

Wi-Fi, USB, Bluetooth, Ethernet, Serial Interface, etc.



Easy to take tissue cassette

U-shaped guide groove design makes it easier for operators to take tissue cassette

Technical parameters

Instrument Name	Tissue Cassette Laser Printer
Instrument Model	TB406I
Printing Technology	Infrared laser printing
Printing Speed	$\leq 3\text{ s}$
Loading Capacity	$\geq 300\text{ pcs}$, 6 slots (optional 450 / 600 pcs)
Output Capacity	$\geq 5\text{ pieces}$
Operation Screen	7-inch capacitive touch screen
Printing Content	Chinese and western characters, numbers, symbols, common one-dimensional codes / two-dimensional codes, logo images
Information Input	Compatible with common hospital interfaces such as LIS / HIS / PIS / PACS
Barcode Scanning	Built-in scanner / external scanner, information is synchronized to the instrument, effectively reducing human input errors
Instrument Size	110 (width) \times 380 length (depth) \times 353 (height) mm
Instrument Weight	13.8 kg



Contact Us
Zenith Biomedical

1265 Barry Ave, PH2, Los Angeles, CA 90025
612-433-7719
info@zenithbiomedical.com