

— Make Pathology Easier

G3061 Slide Laser Printer

Designed for standard slide printing, it uses laser technology to print slide. The instrument is easy to operate, and the printed content is high quality and corrosion-resistant, suitable for long-term storage.



Small • Intelligent • Convenient





High-quality printing

The printed content is clear and can withstand water and xylene. It can be soaked in alcohol and other solutions and is suitable for long-term storage

The template can be edited freely

Free drag and drop editing, arbitrary printing in the slide printing area, support for multiple templates



Dual slots for easy loading

The double slots can hold two different types of slides, loading slides is simple and convenient, improving operating efficiency

Exquisite and compact

The instrument is small in size, with a width of $\leq 14.2\text{cm}$, which greatly saves space




Data transmission

Wi-Fi, USB, Bluetooth, Ethernet, Serial Interface, etc.



Highly compatible slides

Slides with various paint sizes can be printed



Technical parameters

| | |
|---------------------|---|
| Instrument Name | Slide Laser Printer |
| Instrument Model | G306I |
| Printing Technology | Infrared laser printing |
| Print Speed | $\leq 3\text{ s}$ |
| Loading Capacity | $\geq 260\text{ pcs}$ ($\geq 130\text{ pcs} \times 2\text{ slots}$) |
| Output Capacity | 1 pcs, $\geq 20\text{ pcs}$ |
| Operation Screen | 7-inch capacitive touch screen |
| Printing Content | Chinese and western characters, numbers, symbols, common one-dimensional codes / two-dimensional codes, logo images |
| Information Input | Compatible with common hospital interfaces such as LIS / HIS / PIS / PACS |
| Barcode Scanning | Built-in scanner / external scanner, information is synchronized to the instrument, effectively reducing human input errors |
| Instrument Size | 142 (width) \times 298 length (depth) \times 315 (height) mm |
| Instrument Weight | 10.8 kg |



Contact Us
Zenith Biomedical

📍 1265 Barry Ave, PH2, Los Angeles, CA 90025
📞 612-433-7719
✉ info@zenithbiomedical.com